

Glory[®]
Diagnostics



Glory[®] Dengue IgG/IgM Rapid Tests

Glory[®]
Diagnostics



Glory[®]Diagnostics
Manufactured in the Spain

Glory[®] Dengue IgG/IgM Rapit Tests

• Fast

10 Minutes Test Time
Whole Blood Test Protocol

• Accurate

82.4% Sensitivity for IgM in Primary Infection
>99.0% Sensitivity for IgG in Secondary Infection

• Affordable





Innovative Single Use Pack offers logistical relief to users with low-consumption of dengue tests
Economical 10 Devices Regular pack offers economy for users with high consumption




Test Procedure

Dengue IgG/IgM Rapid Test Device


- Add specimen and buffer
- Read results at 10 minutes



IgG & IgM Positive



IgG Positive



IgM Positive



Negative



Invalid

Ordering Information

Catalog No.	Product Name	Kit Size
IDE-402	Glory Dengue IgG/IgM Test Device (WB/S/P)	10 Tests



Authorized Distributor